

## FICHE TECHNIQUE



### **P0142 : PENICILLIN G SODIUM**

CAS number 69-57-8

Molecular weight C<sub>16</sub>H<sub>17</sub>N<sub>2</sub>SO<sub>4</sub>Na = 356.4

Assay 96 %

#### **Storage condition**

- Store protected from light at < 30 °C
- Soluble in water
- Benzylpenicillin sodium

#### **Hazard statements**

- **H317** May cause an allergic skin reaction
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled

#### **Precautionary statements**

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P311** Call a POISON CENTER or doctor/physician
- **P342** If experiencing respiratory symptoms:
- **P352** Wash with soap and water

P0142.0005/ 5 g

P0142.0025/ 25 g

P0142.0100/ 100 g