

## FICHE TECHNIQUE



### **P0159 : PHOSPHINOTRICIN**

CAS number 77182-82-2

Molecular weight  $C_5H_{15}N_2O_4P = 198.2$

#### **Storage condition**

- Store at room temperature
- Soluble in water

#### **Hazard statements**

- **H302** Harmful if swallowed
- **H312** Harmful in contact with skin
- **H332** Harmful if inhaled
- **H360** May damage fertility or the unborn child
- **H373** May cause damage to organs through prolonged or repeated exposure

#### **Precautionary statements**

- **P201** Obtain special instructions before use
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention

P0159.0250/ 250 mg

P0159.1000/ 1 g