

## FICHE TECHNIQUE



### **P0517: POTASSIUM HYDROXIDE**

CAS number 1310-58-3  
Molecular weight KOH = 56.1  
Assay > 85 %

#### **Storage condition**

- Store dry at room temperature
- Soluble in water (1120 g/l / 20 °C)

#### **Hazard statements**

- **H302** Harmful if swallowed
- **H314** Causes severe skin burns and eye damage

#### **Precautionary statements**

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P305** IF IN EYES:
- **P330** Rinse mouth
- **P331** Do NOT induce vomiting
- **P351** Rinse continuously with water for several minutes

P0517.0500/ 500 g

P0517.1000/ 1 kg