

FICHE TECHNIQUE



P0519: POTASSIUM NITRATE

CAS number 7757-79-1
Molecular weight KNO₃ = 101.1
Assay > 99 %

Storage condition

- Store dry at room temperature
- Hygroscopic
- Soluble in water (320 g/l / 20 °C)

Hazard statements

- **H272** May intensify fire; oxidizer

Precautionary statements

- **P210** Keep away from heat/sparks/open flames/hot surfaces

P0519.1000/ 1 kg

P0519.5000/ 5 kg

P0519.9025/ 25 kg