

FICHE TECHNIQUE



P0914 : PICLORAM

CAS number 1918-02-1

Molecular weight $C_6H_3Cl_3N_2O_2 = 241.5$

Storage condition

- Store at room temperature
- Soluble in 1 N NaOH

Hazard statements

- **H319** Causes serious eye irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

P0914.0005/ 5 g

P0914.0010/ 10 g

P0914.0050/ 50 g

RCS 461 800 559 B Montpellier - SIRET 461 800 559 00035 - CODE T.V.A CEE. : FR 62 461 800 559
contact@groupe-meridis.fr - www.groupe-meridis.fr