

Safety Data Sheet

according to Regulation (EU) 2015/830

P0927

Version: 2.0

Date of issue: 28/09/2011 Revision date: 13/11/2017

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Substance

Trade name : Putrescine dihydrochloride

EC-No. : 206-375-9 CAS-No. : 333-93-7 Product code : P0927

Formula : C4H12N2 · 2HCl

Synonyms : 1,4-Diaminobutane / 1,4-Diaminobutanedihydrochloride /

Tetramethylenediaminedihydrochloride

Product group : Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

1.4. Emergency telephone number

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres		

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Serious eye damage/eye irritation, Category 2

Skin corrosion/irritation, Category 2

H315

Specific target organ toxicity — Single exposure, Category 3,

H335

Respiratory tract irritation

Full text of hazard classes and H-statements : see section 16

Safety Data Sheet

according to Regulation (EU) 2015/830

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS07

Signal word (CLP) : Warning

Hazard statements (CLP) : H315 - Causes skin irritation

H319 - Causes serious eye irritation. H335 - May cause respiratory irritation.

Precautionary statements (CLP) : P280 - Wear eye protection, face protection, protective clothing, protective

gloves.

P302+P352 - IF ON SKIN: Wash with plenty of soap and water.

P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name : Putrescine dihydrochloride

CAS-No. : 333-93-7 EC-No. : 206-375-9

Name	Product identifier	%
Putrescine dihydrichloride	(CAS-No.) 333-93-7	> 99
	(EC-No.) 206-375-9	

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Assure fresh air breathing. Get medical advice and attention if you feel

unwell.

First-aid measures after skin contact : Wash skin with plenty of water and soap. Seek medical attention if ill effect

or irritation develops.

First-aid measures after eye contact : Rinse immediately with plenty of water. Seek medical advice.

First-aid measures after ingestion : Rinse mouth thoroughly with water. Seek medical attention if ill effect or

irritation develops.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after skin contact : Redness, pain.

Symptoms/effects after eye contact : Redness, pain. Tears.

13/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: P0927 2/8

Safety Data Sheet

according to Regulation (EU) 2015/830

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

: Alcohol resistant foam. Dry chemical powder. Carbon dioxide (CO2). Water Suitable extinguishing media

spray.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : When exposed to heat, may decompose liberating hazardous gases. COx.

NOx.

5.3. Advice for firefighters

case of fire

Protection during firefighting : Do not enter fire area without proper protective equipment, including

respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Measures in case of dust release : Avoid dust formation.

6.1.2. For emergency responders

No additional information available

6.2. **Environmental precautions**

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Handle in accordance with good industrial hygiene and safety practice.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, well-ventilated area. Store at room temperature. Keep

container closed when not in use.

Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. **Control parameters**

Additional information : Provide local exhaust or general room ventilation

8.2. **Exposure controls**

13/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: P0927 3/8

Safety Data Sheet

according to Regulation (EU) 2015/830

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

Wear suitable protective clothing. Avoid repeated or prolonged skin contact

Respiratory protection:

Wear appropriate mask. Filter type P2 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid

Appearance : Crystals. Powder.

Molecular mass : 161,07 g/mol

Colour : White. Colourless.

Odour : Amine-like.

Odour threshold : No data available pH : No data available Relative evaporation rate : No data available

(butylacetate=1)

Melting point : 280 °C

Freezing point : No data available
Boiling point : No data available
Flash point : No data available

Auto-ignition temperature : ≈ 420

Decomposition temperature : No data available Flammability (solid, gas) : No data available Vapour pressure : No data available Relative vapour density at 20 °C : No data available Relative density : No data available Density : 0,88 g/cm³ (25 °C) Solubility : Water: ≈ 930

Log Pow : -0,84 (25 °C)
Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available

Explosive limits : 0,7 vol %

11,2 vol %

9.2. Other information

No additional information available

13/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: P0927 4/8

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Stable under normal conditions of storage, handling and use.

10.3. Possibility of hazardous reactions

Stable under normal conditions.

10.4. Conditions to avoid

Moisture.

10.5. Incompatible materials

Strong oxidizers.

10.6. Hazardous decomposition products

According to process conditions, hazardous decomposition products may be generated. COx. NOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

Putrescine dihydrochloride (333-93-7)	
LD50 dermal rabbit 825 mg/kg	
LC50 inhalation rat (Dust/Mist - mg/l/4h)	1131 mg/l/4h

Skin corrosion/irritation : Causes skin irritation.

Serious eye damage/irritation : Causes serious eye irritation.

Respiratory or skin sensitisation : Not classified Germ cell mutagenicity : Not classified Carcinogenicity : Not classified

Reproductive toxicity : Not classified

STOT-single exposure : May cause respiratory irritation.

STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

Other information : See actual entry in RTECS for complete information: EJ7280000.

SECTION 12: Ecological information

12.1. Toxicity

Putrescine dihydrochloride (333-93-7)	
LC50 fish 1	730 mg/l Poecilia reticulata (Guppy)

12.2. Persistence and degradability

Putrescine dihydrochloride (333-93-7)	
Biodegradation	>= 89 % (28 d)

Safety Data Sheet

according to Regulation (EU) 2015/830

12.3. Bioaccumulative potential

Putrescine dihydrochloride (333-93-7)	
Log Pow	-0,84 (25 °C)

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Prevent entry to sewers and public waters. Avoid release to the environment

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods

: Dispose of this material and its container at hazardous or special waste collection point. Dispose in a safe manner in accordance with local/national regulations.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA
14.1. UN number		
Not regulated	Not regulated	Not regulated
14.2. UN proper shipping name		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
14.3. Transport hazard class(es)		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
14.4. Packing group		
Not regulated	Not regulated	Not regulated
14.5. Environmental hazards		
Not regulated	Not regulated	Not regulated
No supplementary information available		

14.6. Special precautions for user

Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

13/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: P0927 6/8

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 15: Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

No REACH Annex XVII restrictions

Putrescine dihydrochloride is not on the REACH Candidate List

Putrescine dihydrochloride is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

12th Ordinance Implementing the Federal Immission Control Act -12.BlmSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende

stoffen

: The substance is not listed

SZW-lijst van mutagene stoffen NIET-limitatieve lijst van voor de

voortplanting giftige stoffen -Borstvoeding

: The substance is not listed : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Vruchtbaarheid

: The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Ontwikkeling

: The substance is not listed

Denmark

Recommendations Danish Regulation : Young people below the age of 18 years are not allowed to use the product

Chemical safety assessment 15.2.

No additional information available

SECTION 16: Other information

Indication of changes:

1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves
9.1	Odour	Added	
11.1	LD50 dermal rabbit	Added	
11.1	LC50 inhalation rat (Dust/Mist - mg/l/4h)	Added	
12.1	LC50 fish 1	Added	

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC

Safety Data Sheet

according to Regulation (EU) 2015/830

IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources

: ECHA (European Chemicals Agency). Supplier.

Full text of H- and EUH-statements:

Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Skin Irrit. 2	Skin corrosion/irritation, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H315	Causes skin irritation
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

13/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: P0927 8/8