

FICHE TECHNIQUE



P1353: PHLOROGLUCINOL

CAS number 108-73-6

Molecular weight C₆H₆O₃ = 126.1

Assay 98 %

Storage condition

- Protect from light
- store at room temperature
- Soluble in water

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

P1353.0025/ 25 g

P1353.0100/ 100 g