

FICHE TECHNIQUE

P1362: POLYOXYETHYLENESORBITAN MONOLAURATE

CAS number 9005-64-5

Molecular weight $C_{58}H_{114}O_{26}$ = 1227.7

Fatty acid composition Lauric acid approximately 50 %

Other fatty acids Myristic, Palmitic and Oleic acids

1 l = 1.08 - 1.12 kg

Storage condition

- Store dry at room temperature
- Soluble in water (100 g/l / 25 °C)
- Density: 1.10 kg/l

P1362.0500/ 500 ml

P1362.1000/ 1 l