

Safety Data Sheet

P8004

Version: 3.0

according to Regulation (EU) 2015/830

Supersedes 25/09/2015

Date of issue: 21/01/2011 Revision date: 25/09/2017

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. **Product identifier**

Product form : Substance

Trade name : Pectolyase Y-23

EC-No. : 232-894-5 CAS-No. : 9033-35-6 Product code : P8004

Synonyms : Powdered pectinase preparation (Aspergillus japonicus)

Product group : Raw material

Other means of identification : Enzyme Commission number: 4.2.2.10

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. **Emergency telephone number**

: Supplier contact information: **Emergency number**

> +31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

Classification of the substance or mixture 2.1.

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Respiratory sensitisation, Category 1 H334

Full text of hazard classes and H-statements : see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

according to Regulation (EU) 2015/830

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS08

Signal word (CLP) : Danger

Hazard statements (CLP) : H334 - May cause allergy or asthma symptoms or breathing difficulties if

inhaled.

Precautionary statements (CLP) : P261 - Avoid breathing dust, fume, gas, mist, spray, vapours.

P304+P341 - IF INHALED: if breathing is difficult, remove victim to fresh air

and keep at rest in a position comfortable for breathing.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

 Comments
 : Assay: > 1000 units/g

 Name
 : Pectolyase Y-23

 CAS-No.
 : 9033-35-6

 EC-No.
 : 232-894-5

Name	Product identifier	%
Pectolyase Y-23	(CAS-No.) 9033-35-6	~ 95

Comments : Loss on drying <5%

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Remove victim to fresh air.

First-aid measures after skin contact : Wash off with plenty of water.

First-aid measures after eye contact : Rinse with plenty of water.

First-aid measures after ingestion : Rinse mouth out with water. Seek medical attention if ill effect or irritation

develops.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Water. Foam.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Carbon oxides (CO, CO2). Under fire conditions, hazardous fumes will be

case of fire

present.

5.3. Advice for firefighters

Precautionary measures fire : Wear proper protective equipment.

Firefighting instructions : Prevent fire fighting water from entering the environment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Avoid dust formation. Wear proper protective equipment.

6.1.1. For non-emergency personnel

No additional information available

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

No additional information available

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, cool, well-ventilated area.

Storage temperature : 2 - 8 °C

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

SDS Ref.: P8004

Safety Data Sheet

according to Regulation (EU) 2015/830

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Protective goggles

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Wear appropriate mask. Filter type P3 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powder.
Colour : Brown.

Odour : Characteristic.
Odour threshold : No data available
pH : No data available
Relative evaporation rate : No data available

(butylacetate=1)

: No data available Melting point Freezing point : No data available : No data available **Boiling point** : No data available Flash point Auto-ignition temperature : No data available Decomposition temperature : No data available Flammability (solid, gas) : No data available Vapour pressure : No data available Relative vapour density at 20 °C : No data available Relative density : No data available

Solubility : Water: Readily soluble in water

Log Pow : No data available
Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

No additional information available

10.5. Incompatible materials

No additional information available

10.6. Hazardous decomposition products

No additional information available

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

Skin corrosion/irritation : Not classified Serious eye damage/irritation : Not classified

Respiratory or skin sensitisation : May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Germ cell mutagenicity : Not classified Carcinogenicity : Not classified

Reproductive toxicity : Not classified STOT-single exposure : Not classified STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Prevent entry to sewers and public waters. Avoid release to the

environment.

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.

Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA	
14.1. UN number			
Not regulated	Not regulated	Not regulated	
14.2. UN proper shipping name	14.2. UN proper shipping name		
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.3. Transport hazard class(es)			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.4. Packing group			
Not regulated	Not regulated	Not regulated	
14.5. Environmental hazards			
Not regulated	Not regulated	Not regulated	
No supplementary information available			

14.6. Special precautions for user

. Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

No REACH Annex XVII restrictions

Pectolyase Y-23 is not on the REACH Candidate List

Pectolyase Y-23 is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

6/8

Safety Data Sheet

according to Regulation (EU) 2015/830

VwVwS Annex reference : Water hazard class (WGK) 1, low hazard to waters (Classification according

to VwVwS, Annex 3; ID No. 6949)

12th Ordinance Implementing the Federal Immission Control Act - 12.BlmSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende

stoffen

: Pectolyase Y-23 is listed

: Pectolyase Y-23 is listed: The substance is not listed

SZW-lijst van mutagene stoffen NIET-limitatieve lijst van voor de voortplanting giftige stoffen –

voortplanting giftige stoffen – Borstvoeding

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen –

Ontwikkeling

: The substance is not listed

Denmark

Recommendations Danish Regulation : Young people below the age of 18 years are not allowed to use the product

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

1.1	Other means of identification	Added	
1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Safety Data Sheet

according to Regulation (EU) 2015/830

Data sources : Supplier. REGULATION (EC) No 1272/2008 OF THE EUROPEAN

PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and

amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:

Resp. Sens. 1	Respiratory sensitisation, Category 1
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

 25/09/2017 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: P8004
 8/8