

## FICHE TECHNIQUE

### **P8004 : PECTOLYASE Y-23**

CAS number 9033-35-6

Molecular weight from *Aspergillus japonicus*

Specific activity Approx. 1,000 Maceration units per gram

#### **Storage condition**

- Store dry at 2-8 °C
- Loss on drying < 5%
- Macerating enzyme

#### **Hazard statements**

- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled

#### **Precautionary statements**

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P304** IF INHALED:
- **P341** If breathing is difficult

P8004.0001/ 1 g

P8004.0005/ 5 g