

FICHE TECHNIQUE



Q0250: QUOIRIN & LEPOIVRE MEDIUM

Micro Elements

mg/l

μM

CoCl₂.6H₂O
0.025

CuSO₄.5H₂O
0.025

FeNaEDTA
36.70

H₃BO₃
6.20

KI
0.08

MnSO₄.H₂O
0.76

Na₂MoO₄.2H₂O
0.25

ZnSO₄.7H₂O
8.60

Macro Elements

mg/l

mM

Ca(NO₃)₂.anhydrous

578.92

KH₂PO₄

270.00

KNO₃

1800.00

MgSO₄

175.79

NH₄NO₃

400.00

Druart. P., Sci. Hort., 12, 339-342, (1980).,Quoirin M. and Lepoivre P., Acta Hort, 78, 437, (1977).,Valles. M., Boxus, Ph., Acta Hort., 212, (1987).

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation
- **H412** Harmful to aquatic life with long lasting effects

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P273** Avoid release to the environment
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

Q0250.0001/ 1 l

Q0250.0010/ 10 l

Q0250.0050/ 50 l