

## FICHE TECHNIQUE



### **Q0251: QUOIRIN & LEPOIVRE MEDIUM INCLUDING VITAMINS**

#### Micro Elements

mg/l

$\mu$ M

CoCl<sub>2</sub>.6H<sub>2</sub>O

0.025

CuSO<sub>4</sub>.5H<sub>2</sub>O

0.025

FeNaEDTA

36.70

H<sub>3</sub>BO<sub>3</sub>

6.20

KI

0.08

MnSO<sub>4</sub>.H<sub>2</sub>O

0.76

Na<sub>2</sub>MoO<sub>4</sub>.2H<sub>2</sub>O

0.25

ZnSO<sub>4</sub>.7H<sub>2</sub>O

8.60

## Macro Elements

mg/l

mM

Ca(NO<sub>3</sub>)<sub>2</sub>.anhydrous

578.92

KH<sub>2</sub>PO<sub>4</sub>

270.00

KNO<sub>3</sub>

1800.00

MgSO<sub>4</sub>

175.79

NH<sub>4</sub>NO<sub>3</sub>

400.00

## Vitamins

mg/l

µM

myo-Inositol

100.00

Thiamine HCl

0.40

Total concentration Micro and Macro elements including -vitamins: 3378.40 mg/l  
Druart. P., Sci. Hort., 12, 339-342, (1980).,Quoirin M. and Lepoivre P., Acta Hort, 78, 437,  
(1977).,Valles. M., Boxus, Ph., Acta Hort., 212, (1987).

## Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation
- **H412** Harmful to aquatic life with long lasting effects

## Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P273** Avoid release to the environment
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

Q0251.0001/ 1 l

Q0251.0010/ 10 l

Q0251.0050/ 50 l

Parc Euromédecine -20, rue Robert Koch - BP 27253 34086 MONTPELLIER Cedex 4 - Tel : +33(0)4 67 10 59 00  
RCS 461 800 559 B Montpellier - SIRET 461 800 559 00035 - CODE T.V.A CEE. : FR 62 461 800 559

[contact@groupe-meridis.fr](mailto:contact@groupe-meridis.fr) - [www.groupe-meridis.fr](http://www.groupe-meridis.fr)