

FICHE TECHNIQUE



R0146 : RIFAMPICIN

CAS number 13292-46-1

Molecular weight C₄₃H₅₈N₄O₁₂ = 823.0

Assay > 97 %

Storage condition

- Store dry at 2-8 °C
- Soluble in dilute acid solution

Hazard statements

- **H302** Harmful if swallowed
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation
- **H361** Suspected of damaging fertility or the unborn child
- **H362** May cause harm to breast-fed children

Precautionary statements

- **P202** Do not handle until all safety precautions have been read and understood
- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P263** Avoid contact during pregnancy/while nursing
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention
- **P351** Rinse continuously with water for several minutes

R0146.0001/ 1 g

R0146.0005/ 5 g

