

FICHE TECHNIQUE

S0148: Streptomycin sulphate

CAS number 3810-74-0

Molecular weight ($C_{21}H_{39}N_7O_{12}$) $\cdot 2.3H_2SO_4$ = 1457.4

Assay > 720 IU / mg

Storage condition

Store at 2-8 °C

Slightly soluble in water

Hazard statements

H302 Harmful if swallowed

H360 May damage fertility or the unborn child

Precautionary statements

P201 Obtain special instructions before use

P235 Keep cool

P281 Use personal protective equipment as required

P301 IF SWALLOWED:

P312 Call a POISON CENTER or doctor/physician if you feel unwell

S0148.0050/ 50 g

S0148.0100/ 100 g