

## FICHE TECHNIQUE



### **S0149: SULPHAMETHOXAZOLE**

CAS number 723-46-6

Molecular weight C<sub>10</sub>H<sub>11</sub>N<sub>3</sub>O<sub>3</sub>S = 253.3

#### **Storage condition**

- Store at room temperature
- Soluble in ethanol

#### **Hazard statements**

- **H315** Causes skin irritation
- **H317** May cause an allergic skin reaction
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

#### **Precautionary statements**

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

S0149.0025/ 25 g

S0149.0100/ 100 g