

FICHE TECHNIQUE

S0411

SCHENK&HILDEBRANDT VITAMIN MIXTURE

Vitamins

mg/l

μM

myo-Inositol

1000.0

5549.39

Nicotinic acid

5.0

40.61

Pyridoxine HCl

0.5

2.43

Thiamine HCl

5.0

14.82

Schenk R.U. and Hildebrandt A.C., Medium and techniques for induction and growth of monocotyledonous and dicotyledonous plant cell cultures. Can. J. Bot. 50, 199 (1972).

S0411.0100/ >100 ml

S0411.0250/ >250 ml