

FICHE TECHNIQUE

S0520: SODIUM CHLORIDE

CAS number 7647-14-5

Molecular weight NaCl = 58.4

Assay > 99 %

Bromide (Br) < 0.005 %

Sulphate (SO₄) < 0.02 %

Phosphate (PO₄) < 0.0025 %

Heavy metals < 0.005 %

Storage condition

- Store at room temperature
- Soluble in water (310 g/l / 20 °C)

S0520.1000/ 1 kg

S0520.5000/ 5 kg