

## FICHE TECHNIQUE



### S0523: SODIUM HYDROXIDE

CAS number 1310-73-2  
Molecular weight NaOH = 40.0  
Assay > 98 %

#### Storage condition

- Store dry at room temperature
- Soluble in water (1090 g/l / 20 °C)

#### Hazard statements

- **H290** May be corrosive to metals
- **H314** Causes severe skin burns and eye damage

#### Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P305** IF IN EYES:
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention
- **P330** Rinse mouth
- **P331** Do NOT induce vomiting
- **P351** Rinse continuously with water for several minutes
- **P406** Store in a corrosive resistant/

S0523.0500/ 500 g

S0523.1000/ 1 kg