

FICHE TECHNIQUE



S0524: SODIUM NITRATE

CAS number 7631-99-4

Molecular weight $\text{NaNO}_3 = 85.0$

Assay > 99 %

Storage condition

- Store at room temperature
- Soluble in water (880 g/l / 20 °C)

Hazard statements

- **H272** May intensify fire; oxidizer
- **H302** Harmful if swallowed
- **H319** Causes serious eye irritation

Precautionary statements

- **P210** Keep away from heat/sparks/open flames/hot surfaces
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

S0524.1000/ 1 kg