

FICHE TECHNIQUE



S0536: SILVER NITRATE

CAS number 7761-88-8

Molecular weight AgNO₃ = 169.9

Assay > 99.8 %

Storage condition

- Store at room temperature
- Soluble in water (2150 g/l / 20°C)
- Protect from light

Hazard statements

- **H272** May intensify fire; oxidizer
- **H290** May be corrosive to metals
- **H314** Causes severe skin burns and eye damage
- **H410** Very toxic to aquatic life with long lasting effects

Precautionary statements

- **P260** Do not breathe dust/fume/gas/mist/vapours/spray
- **P273** Avoid release to the environment
- **P301** IF SWALLOWED:
- **P304** IF INHALED:
- **P305** IF IN EYES:
- **P330** Rinse mouth
- **P331** Do NOT induce vomiting
- **P340** Remove victim to fresh air and keep at rest in a position comfortable for breathing
- **P351** Rinse continuously with water for several minutes
- **P370** In case of fire:
- **P378** Use

S0536.0005/ 5 g

S0536.0025/25 g

S0536.0100/ 100 g