

## FICHE TECHNIQUE

### **S1202: SEAPLAQUE™ AGAROSE**

CAS number 9012-36-6  
Gel strength 1.0 % > 200g/cm<sup>2</sup>  
Gelling temperature 26 - 30 C.  
Melting temperature < 65 C.  
Sulphate < 0.10 %  
Elektroendosmosis 0.10  
Moisture < 10 %

#### **Storage condition**

- Store at room temperature
- Soluble in water
- Gel strength: > 200 g/cm<sup>2</sup>

S1202.0100/ 100 g  
S1202.0250/ 250 g