

FICHE TECHNIQUE



S1367: SALICYLIC ACID

CAS number 69-72-7

Molecular weight C₇H₆O₃ = 138.1

Assay > 99 %

Storage condition

- Store dry at room temperature
- Soluble in water (1.8 g/l / 20 °C)

Hazard statements

- **H302** Harmful if swallowed
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention
- **P351** Rinse continuously with water for several minutes

S1367.0100/ 100 g

S1367.0500/ 500 g