

FICHE TECHNIQUE



T0153 : TOBRAMYCINE SULPHATE

CAS number 49842-07-1

Molecular weight (C₁₈H₃₇N₅O₉)₂.5H₂SO₄ = 1425.4

Assay > 634 microg / mg

Storage condition

- Store dry at 2-8 °C
- Soluble in water (50 g/l / 20 °C)

Hazard statements

- **H312** Harmful in contact with skin
- **H332** Harmful if inhaled
- **H360** May damage fertility or the unborn child

Precautionary statements

- **P201** Obtain special instructions before use
- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention

T0153.0001/ 1 g

T0153.0005/ 5 g