

## FICHE TECHNIQUE



### **T0154 : TRIMETHOPRIM**

CAS number 738-70-5

Molecular weight C<sub>14</sub>H<sub>18</sub>N<sub>4</sub>O<sub>3</sub> = 290.3

Assay > 98.5 %

#### **Storage condition**

- Store dry at room temperature
- Soluble in propyleneglycol

#### **Hazard statements**

- **H302** Harmful if swallowed

#### **Precautionary statements**

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell

T0154.0005/ 5 g

T0154.0025/ 25 g