

FICHE TECHNIQUE

T0180 : TICARCILLIN DISODIUM

CAS number 4697-14-7

Molecular weight $C_{15}H_{14}N_2Na_2O_6S_2 = 428.4$

Assay > 95 %

Storage condition

- Store dry at 2-8 °C
- Soluble in water
- hygroscopic

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P311** Call a POISON CENTER or doctor/physician
- **P342** If experiencing respiratory symptoms:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

T0180.0001/ 1 g

T0180.0010/ 10 g