

Draduat identifian

Safety Data Sheet

Date of issue: 14/09/2011

Supersedes 25/09/2015

according to Regulation (EU) 2015/830

T0180

Version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier	
Product form	: Substance
Trade name	: Ticarcillin disodium
CAS-No.	: 4697-14-7
Product code	: T0180
Formula	: C15H14N2Na2O6S2
Synonyms	: Ticarpen

Revision date 05/07/2018

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category	: Professional use
Industrial/Professional use spec	For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. Emergency telephone number

Emergency number

: Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00)

+31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin corrosion/irritation, Category 2	H315
Serious eye damage/eye irritation, Category 2	H319
Respiratory sensitisation, Category 1	H334
Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation	H335

Full text of H statements : see section 16

Adverse physicochemical, human health and environmental effects

Safety Data Sheet

according to Regulation (EU) 2015/830

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)	GHS07 GHS08
Signal word (CLP)	: Danger
Hazard statements (CLP)	 H315 - Causes skin irritation. H319 - Causes serious eye irritation. H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled. H335 - May cause respiratory irritation.
Precautionary statements (CLP)	 P261 - Avoid breathing dust. P280 - Wear protective gloves/protective clothing/eye protection/face protection. P302+P352 - IF ON SKIN: Wash with plenty of water. P342+P311 - If experiencing respiratory symptoms: Call a POISON CENTER/doctor P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name	: Ticarcillin disodium	
CAS-No.	: 4697-14-7	
Name	Product identifier	
Ticarcillin disodium	(CAS-No.) 4697-14-7	

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general	: If you feel unwell, seek medical advice.
First-aid measures after inhalation	: Move to fresh air.
First-aid measures after skin contact	: Wash off with soap and plenty of water
	Seek medical attention if irritation develops.
First-aid measures after eye contact	: Rinse immediately with plenty of water
	Get medical advice/attention.
First-aid measures after ingestion	: Rinse mouth thoroughly with water.

4.2. Most important symptoms and effects, both acute and delayed No additional information available

%

> 95

Safety Data Sheet

according to Regulation (EU) 2015/830

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting m	easures
5.1. Extinguishing media	
Suitable extinguishing media	 Dry chemical powder. Alcohol resistant foam. Carbon dioxide (CO2). Water spray.
5.2. Special hazards arising fro	m the substance or mixture
Hazardous decomposition products in case of fire	 In the second state of the second
5.3. Advice for firefighters Firefighting instructions	: Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Measures in case of dust release : Avoid dust formation. Do not breathe dust.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Do not allow product to spread into the environment.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up

: Sweep up dry powder and dispose properly. Clean contaminated surfaces with an excess of water.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1.	Precautions for safe handling	g
Precau	tions for safe handling	: Handle in accordance with good industrial hygiene and safety practice.
7.2. Conditions for safe storage, including any incompatibilities		
Storage	e conditions	: Store in dry, cool, well-ventilated area

Storage conditions	: Store in dry, cool, well-ventilated area		
	Keep container tightly closed.		
	Moisture sensitive.		
Storage temperature	: 2 - 8 °C		

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Wear approved mask. Filter type P2 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state Appearance Colour	: Solid : Powder. : White.
Odour Odour threshold pH Relative evaporation rate (butylacetate=1)	 No data available No data available No data available No data available
Melting point Freezing point Boiling point Flash point Auto-ignition temperature Decomposition temperature Flammability (solid, gas) Vapour pressure Relative vapour density at 20 °C Relative density Solubility	 No data available So data available No data available No data available So data available No data available No data available
Log Pow Viscosity, kinematic Viscosity, dynamic Explosive properties Oxidising properties Explosive limits	 No data available

Safety Data Sheet

according to Regulation (EU) 2015/830

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions of storage, handling and use.

10.2. Chemical stability

Stable under recommended handling and storage conditions (see section 7).

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

Moisture.

10.5. Incompatible materials

Strong oxidizers.

10.6. Hazardous decomposition products

In the event of fire, may decompose :

- COx
- NOx
- SOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

: Not classified

Ticarcillin disodium (4697-14-7)		
LD50 oral rat	16000 mg/kg	
Skin corrosion/irritation	: Causes skin irritation.	
Serious eye damage/irritation	: Causes serious eye irritation.	
Respiratory or skin sensitisation	: May cause allergy or asthma symptoms or breathing difficulties if inhaled.	
Germ cell mutagenicity	: Not classified	
Carcinogenicity	: Not classified	
Reproductive toxicity	: Not classified	
STOT-single exposure	: May cause respiratory irritation.	
STOT-repeated exposure	: Not classified	
Aspiration hazard	: Not classified	
Other information	: See actual entry in RTECS for complete information: XM9410000.	

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information

: Avoid release to the environment. Prevent entry to sewers and public waters

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods

: Dispose in a safe manner in accordance with local/national regulations.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	ΙΑΤΑ		
14.1. UN number				
Not regulated	Not regulated	Not regulated		
14.2. UN proper shipping name	14.2. UN proper shipping name			
Not regulated	Not regulated	Not regulated		
Not regulated	Not regulated	Not regulated		
14.3. Transport hazard class(es)				
Not regulated	Not regulated	Not regulated		
Not regulated	Not regulated	Not regulated		
14.4. Packing group				
Not regulated	Not regulated	Not regulated		
14.5. Environmental hazards				
Not regulated	Not regulated	Not regulated		
No supplementary information available				

14.6. Special precautions for user

Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture 15.1.1. EU-Regulations

No REACH Annex XVII restrictions Ticarcillin disodium is not on the REACH Candidate List Ticarcillin disodium is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV	: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)
Netherlands	
SZW-lijst van kankerverwekkende stoffen	: The substance is not listed
SZW-lijst van mutagene stoffen	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: The substance is not listed
Denmark	
Recommendations Danish Regulation	: Young people below the age of 18 years are not allowed to use the product

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

	U		
1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et
			cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic

Safety Data Sheet

according to Regulation (EU) 2015/830

REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources

: Supplier. Manufacturer.

Full text of H- and EUH-statements:

Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Resp. Sens. 1	Respiratory sensitisation, Category 1
Skin Irrit. 2	Skin corrosion/irritation, Category 2
Skin Sens. 1	Skin sensitisation, Category 1
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.

SDS Biochemicals version 2018

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product