

## FICHE TECHNIQUE



### **T0190: TICARCILLIN DISODIUM / CLAVULANATE POTASSIUM** Molecular weight MIXTURE

#### **Storage condition**

- Store dry at 2-8 °C
- Soluble in water, Hygroscopic
- CAS No: Ticarcillin 2Na: 4697-14-7 en K Clav.: 61177-45-5

#### **Hazard statements**

- **H315** Causes skin irritation
- **H317** May cause an allergic skin reaction
- **H319** Causes serious eye irritation
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled

#### **Precautionary statements**

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P311** Call a POISON CENTER or doctor/physician
- **P342** If experiencing respiratory symptoms:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

T0190.0002/ 2 g

T0190.0010/ 10 g

T0190.0025/ 25 g