

Safety Data Sheet according to Regulation (EU) 2015/830

T0720

Version: 2.0

Date of issue: 14/09/2011 Revision date 28/02/2018

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Substance
Trade name : L-Tryptophan
EC-No. : 200-795-6
CAS-No. : 73-22-3

REACH registration No : 01-2120006711-71-xxxx

Product code : T0720

Formula : C11H12N2O2

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

1.4. Emergency telephone number

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

SDS Ref.: T0720

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres		

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

Adverse physicochemical, human health and environmental effects

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. **Substances**

Name : L-Tryptophan CAS-No. : 73-22-3 EC-No. : 200-795-6

Name	Product identifier	%
L-Tryptophan	(CAS-No.) 73-22-3	> 98,5
	(EC-No.) 200-795-6	

Full text of H-statements: see section 16

3.2. **Mixtures**

Not applicable

SECTION 4: First aid measures

4.1. **Description of first aid measures**

First-aid measures after inhalation : Assure fresh air breathing

If experiencing respiratory symptoms: Get medical advice/attention.

First-aid measures after skin contact

: Wash skin with mild soap and water.

First-aid measures after eye contact

: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

Obtain medical attention if pain, blinking or redness persists.

First-aid measures after ingestion

: Rinse mouth out with water

Get medical advice/attention if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. **Extinguishing media**

Suitable extinguishing media : Dry chemical, CO2, or water spray or regular foam.

Special hazards arising from the substance or mixture

Hazardous decomposition products in : Under fire conditions, hazardous fumes will be present

case of fire - COx

- NOx.

Safety Data Sheet

according to Regulation (EU) 2015/830

5.3. Advice for firefighters

Firefighting instructions : Prevent fire fighting water from entering the environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including

respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Protective equipment : Wear recommended personal protective equipment.

Measures in case of dust release : Avoid dust formation.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Do not allow to enter drains or water courses. Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Handle in accordance with good industrial hygiene and safety procedures.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, well-ventilated area

Keep container tightly closed Store at ambient temperature Light and moisturer sensitive.

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses

Skin and body protection:

SDS Ref.: T0720

Safety Data Sheet

according to Regulation (EU) 2015/830

Wear suitable protective clothing

Respiratory protection:

Where excessive dust may result, wear approved mask. Type P1 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid

Appearance : Crystalline powder.

Molecular mass : 204,23 g/mol

Colour : White.

Odour : Odourless.

Odour threshold : No data available pH : No data available pH solution : 5,5 - 7 (10 g/l, 20 °C) Relative evaporation rate : No data available

(butylacetate=1)

Melting point : No data available
Freezing point : No data available
Boiling point : No data available
Flash point : No data available
Auto-ignition temperature : No data available

Decomposition temperature : ≈ 290 °C

Flammability (solid, gas) : No data available
Vapour pressure : No data available
Relative vapour density at 20 °C : No data available
Relative density : No data available
Solubility : Water: 10 g/l

Log Pow : -1,05

Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

Sublimation point : 281 - 282 °C (0.4 hPa)

Bulk density : $\approx 400 \text{ kg/m}^3$

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions.

10.2. Chemical stability

Stable under normal conditions of storage, handling and use.

10.3. Possibility of hazardous reactions

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

10.4. Conditions to avoid

- Moisture
- Light (daylight).

10.5. Incompatible materials

- Strong oxidizers.

10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes:

- COx
- NOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

Acute toxicity	. Not classified
L-Tryptophan (73-22-3)	
LD50 oral rat	16 g/kg
Skin corrosion/irritation	: Not classified
Serious eye damage/irritation	: Not classified
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
STOT-single exposure	: Not classified
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified
Other information	: See actual entry in RTECS for complete information: YN6130000.

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

L-Tryptophan (73-22-3)	
Log Pow	-1,05
Bioaccumulative potential	No bioaccumulation.

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

SDS Ref.: T0720

Safety Data Sheet

according to Regulation (EU) 2015/830

12.6. Other adverse effects

Additional information : Avoid release to the environment. Prevent entry to sewers and public waters

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.

SECTION 14: Transport information

In accordance with ADR / IMDG / IATA

ADR	IMDG	IATA	
14.1. UN number			
Not regulated	Not regulated	Not regulated	
14.2. UN proper shipping name			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.3. Transport hazard class(es)			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.4. Packing group			
Not regulated	Not regulated	Not regulated	
14.5. Environmental hazards			
Not regulated	Not regulated	Not regulated	
	No supplementary information available		

14.6. Special precautions for user

. Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

No REACH Annex XVII restrictions

L-Tryptophan is not on the REACH Candidate List

L-Tryptophan is not on the REACH Annex XIV List

15.1.2. National regulations

Germany

Safety Data Sheet

according to Regulation (EU) 2015/830

Reference to AwSV : Water hazard class (WGK) 1, low

hazard to waters (Classification according to VwVwS, Annex 1 or 2;

ID No. 6704)

12th Ordinance Implementing the Federal Immission Control Act -

12.BImSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende stoffen

SZW-lijst van mutagene stoffen

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling

: The substance is not listed : The substance is not listed

: The substance is not listed

: The substance is not listed

: The substance is not listed

Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

1.1	REACH registration No	Added	
1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves
9.2	Sublimation point	Added	
9.2	Bulk density	Added	

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate		
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road		
BCF	Bioconcentration factor		
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008		
DPD	Dangerous Preparations Directive 1999/45/EC		
DSD	Dangerous Substances Directive 67/548/EEC		
IATA	International Air Transport Association		
IMDG	International Maritime Dangerous Goods		
LC50	Median lethal concentration		
LD50	Median lethal dose		
LOAEL	Lowest Observed Adverse Effect Level		
NOAEC	No-Observed Adverse Effect Concentration		
PBT	Persistent Bioaccumulative Toxic		
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006		
SDS	Safety Data Sheet		

Data sources

: Supplier. ECHA (European Chemicals Agency). REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

SDS Biochemicals version 2018

28/02/2018 Duchefa Biochemie B.V. SDS Ref.: T0720 EN (English) 7/8

Safety Data Sheet

according to Regulation (EU) 2015/830

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

 28/02/2018 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: T0720
 8/8