

## FICHE TECHNIQUE



### **T0916 : THIDIAZURON**

CAS number 51707-55-2  
Molecular weight C<sub>9</sub>H<sub>8</sub>N<sub>4</sub>O<sub>S</sub> = 220.2

#### **Storage condition**

- Store at room temperature
- Soluble in DMSO

#### **Hazard statements**

- **H312** Harmful in contact with skin
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H332** Harmful if inhaled
- **H335** May cause respiratory irritation

#### **Precautionary statements**

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

T0916.0250/ 250 mg

T0916.0500/ 500 mg

T0916.1000/ 1 g