

## FICHE TECHNIQUE

### T0929: 2,3,5-TRIIODOBENZOIC ACID

CAS number 88-82-4

Molecular weight  $C_7H_3I_3O_2 = 499.8$

Assay > 95 %

#### Storage condition

- Store dry at 2-8 °C
- Soluble in 1 N NaOH
- Protect from light and moisture

#### Hazard statements

- **H302** Harmful if swallowed

#### Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell

T0929.0005/ 5 g

T0929.0010/ 10 g