

FICHE TECHNIQUE



T1501: TRIS ULTRAPURE

CAS number 77-86-1

Molecular weight $C_4H_{11}NO_3 = 121.1$

A highly purified quality of TRIS with excellent properties for molecular biology and biological buffers. > 98.5 %

Purity, dried substance > 99.9 %

pH (1 M in water) 10.5 - 11.5

Storage condition

- Store dry at room temperature
- Soluble in water (700 g/l / 25 °C)
- Tris(hydroxymethyl)aminomethane

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

T1501.1000/ 1 kg

T1501.5000/ 5 kg