

FICHE TECHNIQUE



T1507: TBE BUFFER

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation
- **H360** May damage fertility or the unborn child

Precautionary statements

- **P202** Do not handle until all safety precautions have been read and understood
- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P281** Use personal protective equipment as required
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention

T1507.0010/ 10 l

T1507.0100/ 100 l