

## Urea

### Safety Data Sheet

according to Regulation (EU) 2015/830 Date of issue: 29/04/2011 Revision date 17/01/2018 Supersedes 02/08/2011

U1363

Version: 2.0

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier	
Product form	: Substance
Trade name	: Urea
EC-No.	: 200-315-5
CAS-No.	: 57-13-6
Product code	: U1363
Formula	: CH4N2O

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

Main use category	: Professional use
Industrial/Professional use spec	: For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

#### **Emergency telephone number** 1.4.

Emergency number

: Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00)

+31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

### **SECTION 2: Hazards identification**

#### Classification of the substance or mixture 2.1.

Classification according to Regulation (EC) No. 1272/2008 [CLP] Not classified

Adverse physicochemical, human health and environmental effects No additional information available

#### 2.2. Label elements

### Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

according to Regulation (EU) 2015/830

### 2.3. Other hazards

No additional information available

## **SECTION 3: Composition/information on ingredients**

### 3.1. Substances

Name	: Urea
CAS-No.	: 57-13-6
EC-No.	: 200-315-5
(	

Name	Product identifier	%
Urea	(CAS-No.) 57-13-6	>= 99
	(EC-No.) 200-315-5	

Full text of H-statements: see section 16

### 3.2. Mixtures

Not applicable

SECTION 4: First aid measures
-------------------------------

### 4.1. Description of first aid measures

: Assure fresh air breathing.
: Wash skin with mild soap and water.
: Rinse cautiously with water for several minutes.
Remove contact lenses, if present and easy to do. Continue rinsing.
: Rinse mouth out with water
Give water to drink.

### 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

## 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures			
<b>5.1. Extinguishing media</b> Suitable extinguishing media Unsuitable extinguishing media	<ul><li>Dry chemical, CO2, or water spray or regular foam.</li><li>Do not use a heavy water stream.</li></ul>		
<b>5.2.</b> Special hazards arising from Hazardous decomposition products in case of fire	the substance or mixture : Thermal decomposition generates : Ammonia (NH3) - COx - NOx.		
<ul><li>5.3. Advice for firefighters</li><li>Firefighting instructions</li><li>Protection during firefighting</li></ul>	<ul><li>Prevent fire fighting water from entering the environment.</li><li>Do not enter fire area without proper protective equipment, including respiratory protection.</li></ul>		

according to Regulation (EU) 2015/830

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

General measures

: Spill should be handled by trained cleaning personnel properly equipped with respiratory and eye protection.

#### 6.1.1. For non-emergency personnel

Protective equipment	: Wear recommended personal protective equipment.
Measures in case of dust release	: Avoid dust formation.

#### 6.1.2. For emergency responders

No additional information available

#### 6.2. Environmental precautions

Avoid release to the environment. Do not allow to enter drains or water courses.

#### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up

: Sweep up dry powder and dispose properly.

### 6.4. Reference to other sections

No additional information available

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

Precautions for safe handling	: Handle in accordance with good industrial hygiene and safety procedures. Ensure good ventilation of the work station.
-------------------------------	--

### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions	: Store in dry, well-ventilated area
	Keep container tightly closed
	Store at room temperature.

### 7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

### **SECTION 8: Exposure controls/personal protection**

#### 8.1. Control parameters

No additional information available

#### 8.2. Exposure controls

#### Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

#### Eye protection:

Safety glasses

#### Skin and body protection:

Long sleeved protective clothing

#### **Respiratory protection:**

Dust production: dust mask with filter type P1

according to Regulation (EU) 2015/830

### SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state Appearance Molecular mass Colour	: Solid : Crystalline powder. : 60,1 g/mol : White.
Odour Odour threshold pH Relative evaporation rate (butylacetate=1)	: ammonia odour. : No data available : < 10 20°C : No data available
Melting point	: 132 - 135 °C
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: ≈ 134 °C
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: No data available
Density	: 1,33 g/cm³ 20°C
Solubility	: Water: ≈ 1000 g/l 20°C
Log Pow	: -1,56
Viscosity, kinematic	: No data available

Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

### 9.2. Other information

No additional information available

### SECTION 10: Stability and reactivity

### 10.1. Reactivity

No additional information available

### 10.2. Chemical stability

Stable under normal conditions.

### **10.3.** Possibility of hazardous reactions

Reacts vigorously with strong oxidizers and acids.

### 10.4. Conditions to avoid

Overheating.

### 10.5. Incompatible materials

No additional information available

U1363

### 10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes: Ammonia (NH3)

- COx

- NOx.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity	: Not classified
Urea (57-13-6)	
LD50 oral rat	14300 mg/kg
Skin corrosion/irritation	: Not classified pH: < 10 20°C
Serious eye damage/irritation	: Not classified pH: < 10 20°C
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
STOT-single exposure	: Not classified
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified

### SECTION 12: Ecological information

### 12.1. Toxicity

Urea (57-13-6)	
LC50 fish 1	6810 mg/l Leuciscus idus (golden orfe)

### 12.2. Persistence and degradability

No additional information available

### 12.3. Bioaccumulative potential

Urea (57-13-6)	
Log Pow	-1,56

### 12.4. Mobility in soil

No additional information available

### 12.5. Results of PBT and vPvB assessment

No additional information available

#### 12.6. Other adverse effects

Additional information

: Avoid release to the environment. Prevent entry to sewers and public waters

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

Waste treatment methods

: Dispose in a safe manner in accordance with local/national regulations.

### **SECTION 14: Transport information**

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA	
14.1. UN number			
Not applicable	Not applicable	Not applicable	
14.2. UN proper shipping name			
Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	
14.3. Transport hazard class(es)			
Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	
14.4. Packing group			
Not applicable	Not applicable	Not applicable	
14.5. Environmental hazards			
Not applicable	Not applicable	Not applicable	
	No supplementary information available		

### 14.6. Special precautions for user

### . Overland transport

Not applicable

### - Transport by sea

Not applicable

#### - Air transport

Not applicable

### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

### **SECTION 15: Regulatory information**

# 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture 15.1.1. EU-Regulations

No REACH Annex XVII restrictions Urea is not on the REACH Candidate List Urea is not on the REACH Annex XIV List

### 15.1.2. National regulations

Germany Reference to AwSV

: Water hazard class (WGK) 1, low hazard to waters (Classification according to VwVwS, Annex 1 or 2; ID No. 118)

## **Urea** Safety Data Sheet

according to Regulation (EU) 2015/830

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV

### Netherlands

SZW-lijst van kankerverwekkende stoffen

SZW-lijst van mutagene stoffen

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting giftige stoffen - Ontwikkeling

### 15.2. Chemical safety assessment

No additional information available

## **SECTION 16: Other information**

### Indication of changes:

		0		
1.4	4	Emergency number	Modified	
2		Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	Obsolete
8.2	2	Hand protection	Modified	Specified material, thickness, et cetera of gloves

### Abbreviations and acronyms:

ATE	Acute Toxicity Estimate	
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road	
BCF	Bioconcentration factor	
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008	
DPD	Dangerous Preparations Directive 1999/45/EC	
DSD	Dangerous Substances Directive 67/548/EEC	
IATA	International Air Transport Association	
IMDG	International Maritime Dangerous Goods	
LC50	Median lethal concentration	
LD50	Median lethal dose	
LOAEL	Lowest Observed Adverse Effect Level	
NOAEC	No-Observed Adverse Effect Concentration	
PBT	Persistent Bioaccumulative Toxic	
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006	
SDS	Safety Data Sheet	

Data sources

: Supplier. ECHA (European Chemicals Agency). REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

### SDS Biochemicals version 2018

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

- : The substance is not listed
- : The substance is not listed