

FICHE TECHNIQUE



V0155: VANCOMYCIN

CAS number 1404-93-9

Molecular weight $C_{66}H_{75}Cl_2N_9O_{24}.HCl=1485.7$

Potency > 1050 IU / mg

Storage condition

- Store dry at 2-8 °C
- Soluble in water (200 g/l)

Hazard statements

- **H317** May cause an allergic skin reaction
- **H319** Causes serious eye irritation

Precautionary statements

- **P235** Keep cool
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P313** Get medical advice/attention
- **P337** If eye irritation persists:
- **P351** Rinse continuously with water for several minutes

V0155.0001/ 1 g

V0155.0005/ 5 g

V0155.0025/ 25 g