

FICHE TECHNIQUE

V0226: VACIN & WENT MEDIUM

Micro Elements

mg/l

µM

Fe₂(C₄H₄O₆)₃

23.13

32.49

MnSO₄.H₂O

5.68

33.61

Macro Elements

mg/l

mM

Ca₃(PO₄)₂

200.00

0.64

KH₂PO₄

250.00

1.84

KNO₃

525.00

5.19

MgSO₄

122.00

1.01

(NH₄)₂SO₄

500.00

3.78

Vacin E.F. and Went E.W., Bot. Gaz. 110, 605 (1949).

V0226.0001/ 1 l

V0226.0010/ 10 l

V0226.0050/ 50 l