

## W1607

### CULTURE TUBES «DE WIT»

Poly carbonate, Gamma radiated  
Height 130 mm  
Diameter middle 27 mm  
Diameter bottom Sterile packed per 75 pieces 10 mm

#### Description :

CULTURE TUBES “DE WIT”,

Polycarbonate, Gamma Radiated Height 130 mm, diameter middle 27 mm, diameter bottom 10 mm.

Culture Tubes “De Wit” are specifically designed for in Vitro Tissue Culture.  
The conical shape of the tubes provides enough space to grow while using a limited quantity of medium.

Culture Tubes “De Wit” are sterile packed per 75 pieces (invoice shows: 10x 75 pcs TOP = W1615 and BOTTOM = W1614)

These culture tubes are not autoclavable! Please beware that Tops & Bottoms are not separate available.

- **W1607 .0750** : 750 pcs
- **W1607 .1500** : 1500 pcs
- **W1607 .2250** : 2250 pcs
- **W1607 .3000** : 3000 pcs
- **W1607 .3750** : 3750 pcs

