

FICHE TECHNIQUE



X1413: X-PHOS P-TOLUIDINE SALT (BCIP P-TOLUIDINE SALT)

CAS number 6578-06-9

Molecular weight $C_8H_6BrClNO_4P.C_7H_9N = 433.6$

Assay > 99 %

Storage condition

- Store at -25 to -15 °C
- Soluble in DMSO and DMF
- Protect from light and moisture

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P411** Store at temperatures not exceeding

X1413.0100/ 100 mg

X1413.1000/ 1 g