

FICHE TECHNIQUE



Z0526: ZINC SULPHATE HEPTAHYDRATE

CAS number 7446-20-0
Molecular weight $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ = 287.6
Assay > 99 %

Storage condition

- Store dry at room temperature
- Soluble in water (960 g/l / 20 °C)

Hazard statements

- **H302** Harmful if swallowed
- **H318** Causes serious eye damage
- **H410** Very toxic to aquatic life with long lasting effects

Precautionary statements

- **P273** Avoid release to the environment
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P305** IF IN EYES:
- **P309** IF exposed or you feel unwell:
- **P310** Immediately call a POISON CENTER or doctor/physician
- **P330** Rinse mouth
- **P331** Do NOT induce vomiting
- **P351** Rinse continuously with water for several minutes

Z0526.0500/ 500 g
Z0526.1000/ 1 kg